

# Imaging of animal communities v1.3

**VIDEO:** [https://www.youtube.com/watch?v=XH\\_dtow-rpg](https://www.youtube.com/watch?v=XH_dtow-rpg)


Imaging has to be done within one month after the extraction to standardize coloration of animals. All animals from each sample are imaged together in 96% ethanol in one Petri dish, 5 photographs per sampling point: nematodes, enchytraeids, mesofauna, macrofauna and earthworms separately (25 photographs per site). All photographs are submitted to the **Central team**, which extracts abundance and individual body masses. This protocol has been developed based on our experience with the scanner EPSON Perfection V600 connected to a PC. Alternative: a scanner which allows for 3200dpi+ and has a CCD image sensor. Petri dish bottom should be flat, thin, and without any scratches (see image).

## Differences between groups:

- Macrofauna and earthworms: plastic Petri dish ø 9 cm, 2400 dpi, white background
- Microarthropods: plastic Petri dish ø 9 cm, 2400 dpi, white background
- Nematoda and enchytraeids: plastic Petri dish ø 3.5 cm, 4800 dpi, black background (the white background panel is removed)

## Key settings:

- Mode: **Professional**
- Export format: **JPEG**
- Document Type: **Reflective**
- Document Source: **Document Table**
- Auto Exposure Type: **Document**
- Image Type: **48-bit color**
- Resolution: **Nematoda and enchytraeids: 4800 dpi; Other: 2400 dpi**
- Unsharp Mask: **on**



## Instructions:

1. Turn on the scanner and the computer
2. Shake the vial with animals and ethanol and pour all animals into a Petri dish
  - a. Ensure to transfer all animals by flushing the vial
  - b. All animals need to be completely covered by ethanol; do not put on the lid
  - c. Use transparent and thin Petri dish without scratches; avoid any bubbles
3. Gently shake the dish back and forth or use tweezers to distribute animals - try to avoid overlap among animals and other objects as much as possible
4. Put the dish into the scanner (in the center). Wait for ~40s so animals sediment.
5. Scanning process
  - a. Create scanning preview and choose scanning frame fitting the Petri dish outer contour (use slightly smaller frame for 9 cm dishes); press "Zoom"
  - b. Press "Scan", choose picture format JPEG (\*.jpg), select the correct folder.
  - c. Name the file according to the **sample code on the label**.
6. Wait until the scanning process is finished *\*\*\*In case multiple scans are taken from the same sample, proceed with step 3.*
7. While the scanning process is running, prepare the next sample.
8. Check the image, put the next dish in the scanner, put the old sample back.